

# Estrogen Receptor alpha / beta Polyclonal Antibody APC Conjugated



Catalog No: #C90092APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Estrogen Receptor alpha / beta Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/595
Conjugates	APC
Target Name	Estrogen Receptor alpha / beta
Other Names	Estrogen Receptor alpha / beta
Accession No.	Uniprot:P03372Gene ID:2099
Uniprot	P03372
GeneID	2099
Calculated MW	66
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only