

Nociceptin receptor Polyclonal Antibody Biotin Conjugated



Catalog No: #C90096Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nociceptin receptor Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1-47/370
Conjugates	Bio
Target Name	Nociceptin receptor
Other Names	Nociceptin receptor
Accession No.	Uniprot:P35377Gene ID:4987
Uniprot	P35377
GeneID	4987
Calculated MW	40
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only