

PDGF-B Polyclonal Antibody Biotin Conjugated

Catalog No: #C90097Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PDGF-B Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	155-210/241
Conjugates	Bio
Target Name	PDGF-B
Other Names	PDGF-B
Accession No.	Uniprot:P01127Gene ID:5155
Uniprot	P01127
GeneID	5155
Calculated MW	27
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only