

# Big endothelin-1 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90098Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Big endothelin-1 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	16-110/212
Conjugates	Bio
Target Name	Big endothelin-1
Other Names	Big endothelin-1
Accession No.	Uniprot:P05305Gene ID:1906
Uniprot	P05305
GeneID	1906
Calculated MW	22
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only