

# alpha smooth muscle Actin Polyclonal Antibody HRP Conjugated



Catalog No: #C90099HRP

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | alpha smooth muscle Actin Polyclonal Antibody HRP Conjugated                    |
| Host Species       | Rabbit  |
| Clonality          | Polyclonal  |
| Purification       | affinity purified by Protein A  |
| Applications       | WB  |
| Species Reactivity | Hu Ms Rat   |
| Immunogen Type     | 301-375/375   |
| Conjugates         | HRP   |
| Target Name        | alpha smooth muscle Actin   |
| Other Names        | alpha smooth muscle Actin   |
| Accession No.      | Uniprot:P62736Gene ID:59  |
| Uniprot            | P62736  |
| GeneID             | 59  |
| Calculated MW      | 42  |
| Formulation        | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                |
| Storage            | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

WB=1:500-2000

Note: This product is for in vitro research use only