

alpha smooth muscle Actin Polyclonal Antibody APC Conjugated



Catalog No: #C90099APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | alpha smooth muscle Actin Polyclonal Antibody APC Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | Flow-Cyt |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 301-375/375 |
| Conjugates | APC |
| Target Name | alpha smooth muscle Actin |
| Other Names | alpha smooth muscle Actin |
| Accession No. | Uniprot:P62736Gene ID:59 |
| Uniprot | P62736 |
| GeneID | 59 |
| Calculated MW | 42 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only