

alpha smooth muscle Actin Polyclonal Antibody PE-Cy5.5 Conjugated



Catalog No: #C90099PE-Cy5.5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	alpha smooth muscle Actin Polyclonal Antibody PE-Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Immunogen Type	301-375/375
Conjugates	PE-Cy5.5
Target Name	alpha smooth muscle Actin
Other Names	alpha smooth muscle Actin
Accession No.	Uniprot:P62736Gene ID:59
Uniprot	P62736
GeneID	59
Calculated MW	42
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only