

## PACAP-27/38 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90100Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PACAP-27/38 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	132-169/176
Conjugates	Bio
Target Name	PACAP-27/38
Other Names	PACAP-27/38
Accession No.	Uniprot:P18509Gene ID:116
Uniprot	P18509
GeneID	116
Calculated MW	19
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only