

GFAP Polyclonal Antibody Biotin Conjugated

Catalog No: #C90105Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GFAP Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	51-150/432
Conjugates	Bio
Target Name	GFAP
Other Names	GFAP
Accession No.	Uniprot:P14136Gene ID:2670
Uniprot	P14136
GeneID	2670
Calculated MW	48
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only