

TSHR (CT) Polyclonal Antibody Biotin Conjugated

Catalog No: #C90108Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TSHR (CT) Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC IF
Species Reactivity	Hu Rat
Immunogen Type	648-740/764
Conjugates	Bio
Target Name	TSHR (CT)
Other Names	TSHR (CT)
Accession No.	Uniprot:P16473Gene ID:7253
Uniprot	P16473
GeneID	7253
Calculated MW	86
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only