

# Hepatitis E Virus ORF3 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90111Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Hepatitis E Virus ORF3 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	HEV
Immunogen Type	22-103/114
Conjugates	Bio
Target Name	Hepatitis E Virus ORF3
Other Names	Hepatitis E Virus ORF3
Accession No.	Uniprot:O90299Gene ID:
Uniprot	O90299
Calculated MW	13.5
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only