

CD95/FAS Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C90112PE-Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD95/FAS Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Ms Rat
Immunogen Type	35-110/327
Conjugates	PE-Cy7
Target Name	CD95/FAS
Other Names	CD95/FAS
Accession No.	Uniprot:P25446Gene ID:246097
Uniprot	P25446
GeneID	246097
Calculated MW	40-50
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only