

HCV Core protein Polyclonal Antibody Biotin Conjugated



Catalog No: #C90115Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HCV Core protein Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	HCV
Immunogen Type	1-50/60
Conjugates	Bio
Target Name	HCV Core protein
Other Names	HCV Core protein
Accession No.	Uniprot:P26664Gene ID:951475
Uniprot	P26664
GeneID	951475
Calculated MW	19/21
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only