

IGF1R Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C90117PE-Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IGF1R Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	251-350/1367
Conjugates	PE-Cy7
Target Name	IGF1R
Other Names	IGF1R
Accession No.	Uniprot:P08069Gene ID:3480
Uniprot	P08069
GeneID	3480
Calculated MW	78/150
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only