

PDGFRA Polyclonal Antibody Biotin Conjugated

Catalog No: #C90120Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PDGFRA Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1021-1089/1089
Conjugates	Bio
Target Name	PDGFRA
Other Names	PDGFRA
Accession No.	Uniprot:P16234Gene ID:5156
Uniprot	P16234
GeneID	5156
Calculated MW	117
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only