

PDGF Receptor beta Polyclonal Antibody Biotin Conjugated

Catalog No: #C90121Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PDGF Receptor beta Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	981-1066/1106
Conjugates	Bio
Target Name	PDGF Receptor beta
Other Names	PDGF Receptor beta
Accession No.	Uniprot:P09619Gene ID:5159
Uniprot	P09619
GeneID	5159
Calculated MW	118
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only