

PLGF Polyclonal Antibody Biotin Conjugated

Catalog No: #C90123Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PLGF Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu Ms
Immunogen Type	131-158/158
Conjugates	Bio
Target Name	PLGF
Other Names	PLGF
Accession No.	Uniprot:P49763Gene ID:18654
Uniprot	P49763
GeneID	18654
Calculated MW	16
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only