

heavy chain cardiac Myosin Polyclonal Antibody PE Conjugated

Catalog No: #C90132PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	heavy chain cardiac Myosin Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1101-1200/1938
Conjugates	PE
Target Name	heavy chain cardiac Myosin
Other Names	heavy chain cardiac Myosin
Accession No.	Uniprot:P12883Gene ID:4624
Uniprot	P12883
GeneID	4624
Calculated MW	213
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only