

Cytomegalovirus pp65 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90135Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytomegalovirus pp65 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC IF
Species Reactivity	Hu HCMVPP65
Immunogen Type	251-350/561
Conjugates	Bio
Target Name	Cytomegalovirus pp65
Other Names	Cytomegalovirus pp65
Accession No.	Uniprot:Q6SW59Gene ID:3077579
Uniprot	Q6SW59
GeneID	3077579
Calculated MW	65
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only