

## Integrin beta 3 Polyclonal Antibody RBITC Conjugated

Catalog No: #C90155RBITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Integrin beta 3 Polyclonal Antibody RBITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms
Immunogen Type	27-120/788
Conjugates	RBITC
Target Name	Integrin beta 3
Other Names	Integrin beta 3
Accession No.	Uniprot:P05106Gene ID:3690
Uniprot	P05106
GeneID	3690
Calculated MW	84
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only