

Protein A Polyclonal Antibody PE Conjugated

Catalog No: #C90162PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein A Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Saureus
Immunogen Type	full length
Conjugates	PE
Target Name	Protein A
Other Names	Protein A
Accession No.	Uniprot:P02976Gene ID:3919448
Uniprot	P02976
GeneID	3919448
Calculated MW	40-60
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only