

beta 2 Microglobulin Polyclonal Antibody Biotin Conjugated

Catalog No: #C90168Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	beta 2 Microglobulin Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F ICC IF
Species Reactivity	Hu
Immunogen Type	21-119/119
Conjugates	Bio
Target Name	beta 2 Microglobulin
Other Names	beta 2 Microglobulin
Accession No.	Uniprot:P61769Gene ID:567
Uniprot	P61769
GeneID	567
Calculated MW	13/27
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only