

MMP14 Polyclonal Antibody Cy3 Conjugated

Catalog No: #C90176Cy3

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | MMP14 Polyclonal Antibody Cy3 Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | ICC IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 271-370/582 |
| Conjugates | Cy3 |
| Target Name | MMP14 |
| Other Names | MMP14 |
| Accession No. | Uniprot:Q9BGL3Gene ID:4323 |
| Uniprot | Q9BGL3 |
| GeneID | 4323 |
| Calculated MW | 54/62 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only