

Streptavidin Polyclonal Antibody Biotin Conjugated

Catalog No: #C90183Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Streptavidin Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-F IF
Species Reactivity	Streptavidin
Conjugates	Bio
Target Name	Streptavidin
Other Names	Streptavidin
Accession No.	Uniprot:Q53532Gene ID:
Uniprot	Q53532
Calculated MW	53
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only