

IL1RAPL1 Polyclonal Antibody FITC Conjugated

Catalog No: #C90187FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL1RAPL1 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	301-400/696
Conjugates	FITC
Target Name	IL1RAPL1
Other Names	IL1RAPL1
Accession No.	Uniprot:Q9NZN1Gene ID:11141
Uniprot	Q9NZN1
GeneID	11141
Calculated MW	78
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only