

GNRH/LHRH Polyclonal Antibody Biotin Conjugated

Catalog No: #C90189Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GNRH/LHRH Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC-P IHC-F IF
Species Reactivity	Hu
Immunogen Type	28-40/92
Conjugates	Bio
Target Name	GNRH/LHRH
Other Names	GNRH/LHRH
Accession No.	Uniprot:P01148Gene ID:2796
Uniprot	P01148
GeneID	2796
Calculated MW	7.9
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only