

Beta casein Polyclonal Antibody FITC Conjugated

Catalog No: #C90191FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Beta casein Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Sheep Go
Immunogen Type	151-222/222
Conjugates	FITC
Target Name	Beta casein
Other Names	Beta casein
Accession No.	Uniprot:P11839Gene ID:443391
Uniprot	P11839
GeneID	443391
Calculated MW	23
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only