

## Ghrelin Polyclonal Antibody RBITC Conjugated

Catalog No: #C90192RBITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ghrelin Polyclonal Antibody RBITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	31-117/117
Conjugates	RBITC
Target Name	Ghrelin
Other Names	Ghrelin
Accession No.	Uniprot:Q9UBU3Gene ID:51738
Uniprot	Q9UBU3
GeneID	51738
Calculated MW	10
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only