

Osteocalcin Polyclonal Antibody PE Conjugated

Catalog No: #C90195PE



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description

|                    |                                                                                 |
|--------------------|---------------------------------------------------------------------------------|
| Product Name       | Osteocalcin Polyclonal Antibody PE Conjugated                                   |
| Host Species       | Rabbit                                                                          |
| Clonality          | Polyclonal                                                                      |
| Purification       | affinity purified by Protein A                                                  |
| Applications       | IF                                                                              |
| Species Reactivity | Hu Ms Rat                                                                       |
| Immunogen Type     | 21-80/100                                                                       |
| Conjugates         | PE                                                                              |
| Target Name        | Osteocalcin                                                                     |
| Other Names        | Osteocalcin                                                                     |
| Accession No.      | Uniprot:P02818Gene ID:632                                                       |
| Uniprot            | P02818                                                                          |
| GeneID             | 632                                                                             |
| Calculated MW      | 11                                                                              |
| Formulation        | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                |
| Storage            | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IF=1:50-200

Note: This product is for in vitro research use only