

SLC1A5 Polyclonal Antibody PE Conjugated

Catalog No: #C90198PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SLC1A5 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu
Immunogen Type	465-515/541
Conjugates	PE
Target Name	SLC1A5
Other Names	SLC1A5
Accession No.	Uniprot:Q15758Gene ID:6510
Uniprot	Q15758
GeneID	6510
Calculated MW	56
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only