

VSIG4 Polyclonal Antibody PE Conjugated

Catalog No: #C90199PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | VSIG4 Polyclonal Antibody PE Conjugated |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | affinity purified by Protein A |
| Applications | IF |
| Species Reactivity | Hu Ms Rat |
| Immunogen Type | 81-160/399 |
| Conjugates | PE |
| Target Name | VSIG4 |
| Other Names | VSIG4 |
| Accession No. | Uniprot:Q9Y279Gene ID:11326 |
| Uniprot | Q9Y279 |
| GeneID | 11326 |
| Calculated MW | 33 |
| Formulation | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

Application Details

IF=1:50-200

Note: This product is for in vitro research use only