

ADRB1 Polyclonal Antibody APC Conjugated

Catalog No: #C90206APC



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ADRB1 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	181-250/477
Conjugates	APC
Target Name	ADRB1
Other Names	ADRB1
Accession No.	Uniprot:P08588Gene ID:153
Uniprot	P08588
GeneID	153
Calculated MW	52
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only