

## CD45 Polyclonal Antibody APC Conjugated

Catalog No: #C90213APC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD45 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Ms Rat
Immunogen Type	1210-1304/1304
Conjugates	APC
Target Name	CD45
Other Names	CD45
Accession No.	Uniprot:P08575Gene ID:5788
Uniprot	P08575
GeneID	5788
Calculated MW	143
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only