

CD24 Polyclonal Antibody FITC Conjugated

Catalog No: #C90217FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD24 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rat
Immunogen Type	60-80/80
Conjugates	FITC
Target Name	CD24
Other Names	CD24
Accession No.	Uniprot:P25063Gene ID:100133941
Uniprot	P25063
GeneID	100133941
Calculated MW	8
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200

Note: This product is for in vitro research use only