

## IL-18 Polyclonal Antibody FITC Conjugated

Catalog No: #C90218FITC

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IL-18 Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Ms Rat
Immunogen Type	51-120/192
Conjugates	FITC
Target Name	IL-18
Other Names	IL-18
Accession No.	Uniprot:P70380Gene ID:16173
Uniprot	P70380
GeneID	16173
Calculated MW	17
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

IF=1:50-200

Note: This product is for in vitro research use only