

PPAR gamma Polyclonal Antibody FITC Conjugated

Catalog No: #C90219FITC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PPAR gamma Polyclonal Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/505
Conjugates	FITC
Target Name	PPAR gamma
Other Names	PPAR gamma
Accession No.	Uniprot:P37231Gene ID:5468
Uniprot	P37231
GeneID	5468
Calculated MW	57
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test IF=1:50-200

Note: This product is for in vitro research use only