

PPAR gamma Polyclonal Antibody Cy5 Conjugated

Catalog No: #C90219Cy5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PPAR gamma Polyclonal Antibody Cy5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	101-200/505
Conjugates	Cy5
Target Name	PPAR gamma
Other Names	PPAR gamma
Accession No.	Uniprot:P37231Gene ID:5468
Uniprot	P37231
GeneID	5468
Calculated MW	57
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only