

Chromogranin A Polyclonal Antibody PE Conjugated

Catalog No: #C90223PE



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Chromogranin A Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	421-470/457
Conjugates	PE
Target Name	Chromogranin A
Other Names	Chromogranin A
Accession No.	Uniprot:P10645Gene ID:1113
Uniprot	P10645
GeneID	1113
Calculated MW	48
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only