

IL-13 Polyclonal Antibody Bio Conjugated

Catalog No: #C90231Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-13 Polyclonal Antibody Bio Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F IF
Species Reactivity	Hu Ms Rat
Conjugates	Bio
Target Name	IL-13
Other Names	IL-13
Accession No.	Uniprot:P20109Gene ID:16163
Uniprot	P20109
GeneID	16163
Calculated MW	12
Concentration	1mg/ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only