

CD62p Polyclonal Antibody Biotin Conjugated

Catalog No: #C90232Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD62p Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F IF
Species Reactivity	Hu Ms Rat
Immunogen Type	701-768/768
Conjugates	Bio
Target Name	CD62p
Other Names	CD62p
Accession No.	Uniprot:Q01102Gene ID:6403
Uniprot	Q01102
GeneID	6403
Calculated MW	88
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only