

VEGFR2 Polyclonal Antibody PE-Cy5.5 Conjugated

Catalog No: #C90236PE-Cy5.5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VEGFR2 Polyclonal Antibody PE-Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	1001-1100/1356
Conjugates	PE-Cy5.5
Target Name	VEGFR2
Other Names	VEGFR2
Accession No.	Uniprot:P35968Gene ID:3791
Uniprot	P35968
GeneID	3791
Calculated MW	147
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only