

Smad4 Polyclonal Antibody PE-Cy5.5 Conjugated

Catalog No: #C90247PE-Cy5.5

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Smad4 Polyclonal Antibody PE-Cy5.5 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IF
Species Reactivity	Hu Ms Rat
Immunogen Type	31-120/545
Conjugates	PE-Cy5.5
Target Name	Smad4
Other Names	Smad4
Accession No.	Uniprot:Q13485Gene ID:4089
Uniprot	Q13485
GeneID	4089
Calculated MW	60
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Note: This product is for in vitro research use only