

## IL-2 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90251Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IL-2 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC-P IHC-F IF
Species Reactivity	Hu
Immunogen Type	21-110/153
Conjugates	Bio
Target Name	IL-2
Other Names	IL-2
Accession No.	Uniprot:P60568Gene ID:3558
Uniprot	P60568
GeneID	3558
Calculated MW	15
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only