

## ICAM1/CD54 Polyclonal Antibody PE Conjugated

Catalog No: #C90252PE

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ICAM1/CD54 Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	201-300/537
Conjugates	PE
Target Name	ICAM1/CD54
Other Names	ICAM1/CD54
Accession No.	Uniprot:P13597Gene ID:3383
Uniprot	P13597
GeneID	3383
Calculated MW	56
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only