

CD272 Polyclonal Antibody PE-Cy7 Conjugated

Catalog No: #C90256PE-Cy7

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD272 Polyclonal Antibody PE-Cy7 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	221-289/289
Conjugates	PE-Cy7
Target Name	CD272
Other Names	CD272
Accession No.	Uniprot:Q7Z6A9Gene ID:151888
Uniprot	Q7Z6A9
GeneID	151888
Calculated MW	28
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only