

Aquaporin 4 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90257Bio

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aquaporin 4 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB IHC P IHC-F ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	271-323/323
Conjugates	Bio
Target Name	Aquaporin 4
Other Names	Aquaporin 4
Accession No.	Uniprot:P55088Gene ID:361
Uniprot	P55088
GeneID	361
Calculated MW	36
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:50-200 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only