

## Integrin alpha 4 Polyclonal Antibody Biotin Conjugated

Catalog No: #C90261Bio

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Integrin alpha 4 Polyclonal Antibody Biotin Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	WB
Species Reactivity	Hu Ms Rat
Immunogen Type	551-650/1032
Conjugates	Bio
Target Name	Integrin alpha 4
Other Names	Integrin alpha 4
Accession No.	Uniprot:P13612Gene ID:3676
Uniprot	P13612
GeneID	3676
Calculated MW	111
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

WB=1:50-200

Note: This product is for in vitro research use only