

C-jun Polyclonal Antibody PE Conjugated

Catalog No: #C90273PE

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-jun Polyclonal Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rat
Immunogen Type	31-331/331
Conjugates	PE
Target Name	C-jun
Other Names	C-jun
Accession No.	Uniprot:P05412Gene ID:3725
Uniprot	P05412
GeneID	3725
Calculated MW	43
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only