

FGFR2 Polyclonal Antibody APC Conjugated

Catalog No: #C90275APC

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FGFR2 Polyclonal Antibody APC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rat
Immunogen Type	21-120/821
Conjugates	APC
Target Name	FGFR2
Other Names	FGFR2
Accession No.	Uniprot:P21802Gene ID:2263
Uniprot	P21802
GeneID	2263
Calculated MW	89
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test ICC=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only