

## CD147 Polyclonal Antibody PE-Cy5 Conjugated

Catalog No: #C90279PE-Cy5

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |                                                                                 |
|--------------------|---------------------------------------------------------------------------------|
| Product Name       | CD147 Polyclonal Antibody PE-Cy5 Conjugated                                     |
| Host Species       | Rabbit                                                                          |
| Clonality          | Polyclonal                                                                      |
| Purification       | affinity purified by Protein A                                                  |
| Applications       | Flow-Cyt IF                                                                     |
| Species Reactivity | Hu Ms Rat                                                                       |
| Immunogen Type     | 301-385/385                                                                     |
| Conjugates         | PE-Cy5                                                                          |
| Target Name        | CD147                                                                           |
| Other Names        | CD147                                                                           |
| Accession No.      | Uniprot:P35613Gene ID:682                                                       |
| Uniprot            | P35613                                                                          |
| GeneID             | 682                                                                             |
| Calculated MW      | 40                                                                              |
| Formulation        | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                |
| Storage            | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Note: This product is for in vitro research use only